

DEPArray™: Progettazione, prototipazione e fabbricazione di un sistema per l'individuazione, la manipolazione ed il sorting di singole cellule

Ing. Gianni Medoro (Silicon Biosystems)
Dott. Davide Ramanand (Skytechnology)

NI Medical Innovation Summit Modena 8th June 2010

Silicon Biosystems Spa – Bologna ITALY www.siliconbiosystems.com

Agenda

- Company Introduction
- ➤ DEPArrayTM Technology
 - Cell Manipulation by Dielectrophoresis
 - Uniqueness of Performances
- ➤ DEPArrayTM Applications
 - Oncology: Personalized Therapy
 - Prenatal Diagnosis: Non Invasive Approach
- Partnership Network
 - Synergy with National Instruments and Sky Technology



10 Years from Technology Concept to Mature Product with Applications

- 1999 Company founded by 4 scientists at Bologna University Core patent on Moving DEP Cages filed
- 2000-2005 DEPArray ™ technology Proof-of-concept developed under EU, Government and Regional funding (> 5M€ total)
 - 2006 1st fund raising round (1.5M€)
 - DEPArray™ technology and NIPD application development
 - 2007 2nd fund raising round (5.5M€)

 Direct & Indirect Gov funding: (~ 2.5 M€)
 - 2008 3rd fund raising round (3.5M€)
 Oncology application development started
 - 2009 CE Mark for DEPArray™
 Start of marketing activities
 - 2010 4th fund raising round (5M€)

 Head count ~24 FTE, 27 patents filed


 NIPD and Oncology Applications close to Pilot Trial

DEPArray™ Cell Sorting


Silicon Biosystems Patented Technology: Moving DEP Cages

Cell trapping by DEP cages

Non-uniform electric field generated by the chip
electrodes (cross section)
cage-move




Sorting of Fluorescent Cells with DEPArray™





Moving DEP Cages Enables Outstanding Performance in Single-Cell Sorting



- 1. Inject, trap and image all cells
- 2. Move all cells of interest into Parking chamber
- 3. Move separately to Recovery chamber and flush

Features

- Multiparametric image-based sorting
- ➤ Single cell resolution
- Small cell loads
- > Gentle on cells
- No a priori thresholds → choose best cells
- Sorting based on slow kinetics (10-100s minutes) possible




Individual Cell Sorting with DEPArray™ A300K



Click on image to start movie



DEPArray™ is the *Last Mile* for Automatic Rare-Cells Sorting in High-Value Applications




Targeting Unmet Medical Needs with DEPArray™

- Oncology: Circulating Tumor Cells Isolation
 - From prognosis (based on enumeration of CTC) to personalized therapy (based on isolation of CTC)
- > Prenatal Diagnosis: Circulating Fetal Cells Isolation
 - From Amniocentesis and Villocentesis to Non-Invasive Prenatal Diagnosis (based on Circulating Fetal Cells)
- ➤ Life-Science Research
 - Stem-cells
 - Immunology
 - Single-cell biology



Cancer Personalized Therapy with DEPArrayTM

Rare, but Important: The Role of CTC in the Creation of Metastasis



The bone marrow niche: habitat to hematopoietic and mesenchymal stem cells, and unwitting host to molecular parasites

Y Shiozawa, A M Havens, K J Pienta and R S Taichman




Personalized Oncology Treatment: genetic detection to identify addressable patients

- Silicon Biosystems tests for genetic profiling enables to determine effectiveness and allows targeting of treatment
 - Cell growth inhibitor treatment is fastest growing class of cancer drugs from \$9.5bn in 2005 to \$22bn in 2010 at a CAGR of 18%
 - Cell growth inhibitor treatment is expensive (monthly cost ranging from \$3,000 to \$10,000), with significant side-effects

Annual courses of treatment EU, US, Japan by type of cancer

Colorectal and Lung represent ca 1m courses of treatment requiring monitoring each year





CTC <u>genetic profiling</u>: dramatic improvement in patient life



A market paradigm switch: from prognosis to personalized therapy

CTC molecular analysis allows selection of most effective drugs


- Longer & Better life for patient
- Cost saving for healthcare system




DEPArray™ Enables the Isolation of Circulating Tumor Cells and their Molecular Analysis

➤ Genome-wide detection of chromosomal

unbalances




CTCs (CK-green/DAPI-blue) isolated with DEPArray™ from a metastatic breast cancer patient



Detection of point mutations (e.g. KRAS)

CTCs isolated with DEPArray™ from a metastatic colorectal cancer patient




e.g. KRAS mutation testing is required to predict efficacy of treatment with Ceutximab or Panitumumab



Non Invasive Prenatal Diagnosis with DEPArray™


Amniocentesis and chorionic villus sampling: high-risk of miscarriages and unpleasant procedures

Amniocentesis



Risk of fetal loss 0.5-1.5%

Chorionic Villus Sampling




Mother's age	Risk	Need
20 years	1 in 1,600	Low
25 years	1 in 1,300	Medium
30 years	1 in 1,000	
35 years	1 in 365	
40 years	1 in 90	High
45 years	1 in 30	



DEPArray™ Enables the Isolation of Circulating Fetal Cells and their Molecular Analysis

> Detection and Sorting of fetal cells from


maternal blood



Example of Fetal Cell isolated with DEPArray[™] from peripheral blood of a pregnant woman:


nRBC (γ-Globin-red/DAPI-blue)

➤ Genetic Analysis





Sizeable and highly predictable market for NIPD



2.0m high-risk pregnancies in EU, US and Japan, but only 24% adoption of amniocentesis test which has a risk of miscarriage



Partnership Network


DEPArray™ System




Equipment



Lab-on-a-Chip



Software





Network of Competences

Mixed Analog-Digital I/O

Chip Manufacturing Partner

Image Processing

Equipment
Manufacturing
Partner



Software &
Firmware IVD
Support

Heat Management

Power RF Design



Network of Competences



Chip Manufacturing Partner



Equipment
Manufacturing
Partner






Heat Management

Power RF Design



Partnership goals





Needs



- Flexibility
- Time-to-market
- Proven success (on specific items)
- Qualified support
- Future development support
- Robustness
- Maintainbility



Why NI - 1



- Worldwide player
 - Products wide choice ('Time-to-market', 'Support')
 - Ease of integration (HW with SW) ('Time-to-Market', 'Flexibility')
 - Local & central branches ('Support')
- Corporate support
 - R&D upon agreement ('Support')
 - Product Development ('Future development')
- Alliance Partner Network
 - Services Firm ('Support', 'Robustness')
 - Integrators ('Support', 'Robustness')
- Large Libraries Collection ('Time-to-market', 'Support', 'Robustness')
- Long-Term Support ('Maintenability' (& Evolution as well))
- LabVIEW real-time ('Robustness')



Why NI - 2



- Try-'n'-buy operation ('Time-to-Market', 'Support' (NI Italy+corporate))
- Openness
 - Integration with Silicon Biosystems proprietary service infrastructure ('Time-to-Market', 'Robustness')
 - Integration with third-parties instrumentation drivers ('Time-to-Market')
 - Integration with third-parties imaging algorithmics ('Time-to-Market', 'Robustness')
- Wide COTS offering
 - Cost optimization path support ((near!) 'Future development') (FlexRIO?)
 - Runtime changes can be easily managed ('Flexibility', 'Time-to-market')
- Certified/Supported COTS
 - Widely used COTS as a facilitating element in regulated markets certification achievement ('Support', 'Robustness')

In our experience, NI is a dynamic reality putting much attention on market evolutions, interested in supporting strategic partnerships.



Conclusions



- Technically: an outstanding exercise!
 - The machine involves knowledge and requirements necessitating many different aspects and corresponding knowledge
- Ethically: biomedical devices are of undiscussed ethic value
 - Thinking in terms of machine intended use is a great motivational point
- But overall: great team work!
 - All involved parties have really pushed far beyond 'the usual'
 - Integration issues couldn't be surpassed without all actors strong involvment

...So, looking forward to next meeting to share with you future evolutions!



www.siliconbiosystems.com

www.ni.com/it

www.skytechnology.it